

## Product datasheet for **AM20802PU-N**

### **SERPING1 Mouse Monoclonal Antibody [Clone ID: KT28]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	KT28
Applications:	ELISA
Recommended Dilution:	ELISA: 1 µg/ml.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Purified C1 inhibitor from human blood
Specificity:	This antibody reacts to SERPING1.
Formulation:	PBS containing 0.09% sodium azide as a preservative State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Protein A chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	serpin family G member 1
Database Link:	<a href="#">Entrez Gene 710 Human P05155</a>
Synonyms:	C1IN, C1NH, C1 Inh, C1Inh
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Complement and coagulation cascades



[View online »](#)