

## Product datasheet for **AM20803PU-N**

### Factor I (CFI) Mouse Monoclonal Antibody [Clone ID: KT22]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	KT22
Applications:	ELISA
Recommended Dilution:	ELISA: 1 µg/ml.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Purified Factor I from human blood
Specificity:	This antibody reacts to Complement factor I.
Formulation:	PBS containing 0.09% sodium azide as a preservative State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Protein A chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	complement factor I
Database Link:	<a href="#">Entrez Gene 3426 Human P05156</a>
Synonyms:	CFI, IF, KAF, AHUS3, C3BINA, C3b-INA
Protein Families:	Druggable Genome, Protease, Secreted Protein
Protein Pathways:	Complement and coagulation cascades



[View online »](#)