

## Product datasheet for **AM20805PU-S**

### Complement factor B (CFB) Mouse Monoclonal Antibody [Clone ID: KT21]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	KT21
Applications:	ELISA
Recommended Dilution:	Indirect ELISA: 1 µg/ml. Sandwich ELISA: 5 µg/ml.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Purified Factor B from human blood
Specificity:	This antibody reacts to Complement factor B.
Formulation:	PBS containing 0.09% sodium azide as a preservative State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Protein A chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	complement factor B
Database Link:	<a href="#">Entrez Gene 629 Human P00751</a>
Synonyms:	Properdin factor B, C3/C5 convertase, PBF2, CFB, BF, BFD
Protein Families:	Druggable Genome, Protease, Secreted Protein
Protein Pathways:	Complement and coagulation cascades



[View online »](#)