

Product datasheet for **AM20806PU-L**

Properdin (CFP) Mouse Monoclonal Antibody [Clone ID: KT20]

Product data:

Product Type:	Primary Antibodies
Clone Name:	KT20
Applications:	ELISA
Recommended Dilution:	Indirect ELISA: 1 µg/ml.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human properdin purified from blood
Specificity:	This antibody reacts to Properdin.
Formulation:	PBS containing 0.09% sodium azide as a preservative State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Protein A chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	complement factor properdin
Database Link:	Entrez Gene 5199 Human P27918
Synonyms:	CFP, PFC, Complement factor P
Protein Families:	Secreted Protein



[View online »](#)