

Product datasheet for **AM20807PU-S**

CD35 (CR1) Mouse Monoclonal Antibody [Clone ID: KT23]

Product data:

Product Type:	Primary Antibodies
Clone Name:	KT23
Applications:	FC
Recommended Dilution:	Flow cytometry: 10 µg/ml.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Recombinant soluble human CR1
Specificity:	This antibody reacts to CD35.
Formulation:	PBS containing 0.09% sodium azide as a preservative State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Protein A chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	complement component 3b/4b receptor 1 (Knops blood group)
Database Link:	Entrez Gene 1378 Human P17927
Synonyms:	Complement receptor type 1, CR1, C3b/C4b receptor
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Complement and coagulation cascades, Hematopoietic cell lineage



[View online »](#)