

## Product datasheet for **AM20808PU-S**

### Human IgA (Fc specific) (IgA1 + IgA2) Mouse Monoclonal Antibody [Clone ID: KT13]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgA (Fc specific) (IgA1 + IgA2) Mouse Monoclonal Antibody [Clone ID: KT13]
Clone Name:	KT13
Applications:	ELISA
Recommended Dilution:	Indirect ELISA: 1 µg/ml.
Reactivity:	Human
Host:	Mouse
Immunogen:	Recombinant Human IgA2
Isotype:	IgG2a
Formulation:	PBS containing 0.09% Sodium Azide as a preservative State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Protein A Chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.



[View online »](#)