

Product datasheet for **AM20812PU-S**

IL2 Mouse Monoclonal Antibody [Clone ID: KT7]

Product data:

Product Type:	Primary Antibodies
Clone Name:	KT7
Applications:	ELISA
Recommended Dilution:	ELISA: 1 µg/ml.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Recombinant human IL-2
Specificity:	This antibody reacts to Interleukin-2.
Formulation:	PBS containing 0.09% sodium azide as a preservative State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Protein A chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	interleukin 2
Database Link:	Entrez Gene 3558 Human P17108
Synonyms:	IL-2, TCGF



[View online »](#)