

## Product datasheet for **AM20813PU-L**

### IL11 Mouse Monoclonal Antibody [Clone ID: KT8]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	KT8
Applications:	ELISA
Recommended Dilution:	ELISA: 1 µg/ml. Western blotting: 1 µg/ml. Sandwich ELISA: 5 µg/ml for coating plates.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Recombinant human IL-11
Specificity:	This antibody reacts to Interleukin-11.
Formulation:	PBS containing 0.09% sodium azide as a preservative State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Protein A chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	interleukin 11
Database Link:	<a href="#">Entrez Gene 3589 Human P20809</a>
Synonyms:	IL-11, Adipogenesis inhibitory factor, AGIF
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)

**Protein Pathways:** Cytokine-cytokine receptor interaction, Hematopoietic cell lineage, Jak-STAT signaling pathway