

Product datasheet for **AM20814PU-S**

IL1RA (IL1RN) Mouse Monoclonal Antibody [Clone ID: KT6]

Product data:

Product Type:	Primary Antibodies
Clone Name:	KT6
Applications:	ELISA
Recommended Dilution:	Indirect ELISA: 1 µg/ml as primary antibody.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Recombinant human IL-1ra.
Specificity:	This antibody reacts to IL1RN.
Formulation:	PBS containing 0.09% sodium azide as a preservative State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Protein A chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	interleukin 1 receptor antagonist
Database Link:	Entrez Gene 3557 Human P18510
Synonyms:	IL1F3, IL-1ra, IL-1RN, IRAP, IL1 inhibitor, ICIL-1RA
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)