

Product datasheet for **AM20816PU-S**

Endostatin (COL18A1) Mouse Monoclonal Antibody [Clone ID: KT3]

Product data:

Product Type:	Primary Antibodies
Clone Name:	KT3
Applications:	ELISA
Recommended Dilution:	ELISA: 1 µg/ml.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Recombinant endostatin
Specificity:	This antibody reacts to Collagen type XVIII alpha 1 chain.
Formulation:	PBS containing 0.09% sodium azide as a preservative State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Protein A chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	collagen type XVIII alpha 1 chain
Database Link:	Entrez Gene 80781 Human P39060
Synonyms:	COL18A1, Endostatin, KNO1



[View online »](#)