

Product datasheet for **AM20819PU-N**

Clasp1 (N-term) Rat Monoclonal Antibody [Clone ID: KT67]

Product data:

Product Type:	Primary Antibodies
Clone Name:	KT67
Applications:	IF, WB
Recommended Dilution:	Western blot: 1 µg/ml. Immunofluorescent staining: 10 µg/ml.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	alpha-CLASP1-N-terminus and GST fusion protein
Specificity:	This antibody reacts to CLASP1 (N-term).
Formulation:	PBS containing 0.09% sodium azide as a preservative State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Protein G chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	CLIP associating protein 1
Database Link:	Entrez Gene 76707 Mouse Q80TV8
Synonyms:	CLIP-associating protein 1, KIAA0622, MAST1, hOrbit1



[View online »](#)