

Product datasheet for **AM20820PU-S**

Rab11b Rat Monoclonal Antibody [Clone ID: KT80]

Product data:

Product Type:	Primary Antibodies
Clone Name:	KT80
Applications:	IF, WB
Recommended Dilution:	Western blotting: 1 µg/ml. Immunofluorescent staining: 10 µg/ml.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Mouse RAB11B and GST fusion protein
Specificity:	This antibody reacts to RAB11B.
Formulation:	PBS containing 0.09% sodium azide as a preservative State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Protein G chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	RAB11B, member RAS oncogene family
Database Link:	Entrez Gene 19326 Mouse P46638
Synonyms:	Rab-11B, YPT3



[View online »](#)