

## Product datasheet for **AM20822PU-L**

### **Rab6b Rat Monoclonal Antibody [Clone ID: KT79]**

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	KT79
Applications:	IF, WB
Recommended Dilution:	Western blotting: 0.1 - 1 µg/ml. Immunofluorescent staining: 1 - 10 µg/ml.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Mouse RAB6B
Specificity:	This antibody reacts to RAB6B.
Formulation:	PBS containing 0.09% sodium azide as a preservative State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Protein G chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	RAB6B, member RAS oncogene family
Database Link:	<a href="#">Entrez Gene 270192 Mouse P61294</a>
Synonyms:	Rab-6B



[View online »](#)