

## Product datasheet for **AM20824PU-L**

### **Chtop (C-term) Rat Monoclonal Antibody [Clone ID: KT64]**

#### **Product data:**

|                       |                                                                                                                                                     |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:         | Primary Antibodies                                                                                                                                  |
| Clone Name:           | KT64                                                                                                                                                |
| Applications:         | IF, IHC, IP, WB                                                                                                                                     |
| Recommended Dilution: | Western Blot: 1/1000.<br>Immunofluorescent staining: 1/100.<br>Immunohistochemistry on frozen sections: 1/100.<br>Immunoprecipitation: 1 - 2 µg/IP. |
| Reactivity:           | Mouse                                                                                                                                               |
| Host:                 | Rat                                                                                                                                                 |
| Isotype:              | IgG2a                                                                                                                                               |
| Clonality:            | Monoclonal                                                                                                                                          |
| Immunogen:            | Mouse Friend of Prmt1 C-terminus and GST fusion protein                                                                                             |
| Specificity:          | This antibody reacts to C1orf77 (C-term).                                                                                                           |
| Formulation:          | PBS containing 0.09% sodium azide as a preservative<br>State: Aff - Purified<br>State: Liquid purified Ig fraction                                  |
| Concentration:        | lot specific                                                                                                                                        |
| Purification:         | Protein G chromatography                                                                                                                            |
| Conjugation:          | Unconjugated                                                                                                                                        |
| Storage:              | Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C to -80°C for longer.<br>Avoid repeated freezing and thawing.          |
| Stability:            | Shelf life: one year from despatch.                                                                                                                 |
| Gene Name:            | chromatin target of PRMT1                                                                                                                           |
| Database Link:        | <a href="#">Entrez Gene 66511 Mouse Q9CY57</a>                                                                                                      |
| Synonyms:             | HT031, PP7704                                                                                                                                       |



[View online »](#)