

## Product datasheet for **AM20824PU-N**

### Chtop (C-term) Rat Monoclonal Antibody [Clone ID: KT64]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	KT64
Applications:	IF, IHC, IP, WB
Recommended Dilution:	Western Blot: 1/1000. Immunofluorescent staining: 1/100. Immunohistochemistry on frozen sections: 1/100. Immunoprecipitation: 1 - 2 µg/IP.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Mouse Friend of Prmt1 C-terminus and GST fusion protein
Specificity:	This antibody reacts to C1orf77 (C-term).
Formulation:	PBS containing 0.09% sodium azide as a preservative State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Protein G chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	chromatin target of PRMT1
Database Link:	<a href="#">Entrez Gene 66511 Mouse Q9CY57</a>
Synonyms:	HT031, PP7704



[View online »](#)