

## Product datasheet for **AM20825PU-N**

### APC / DP2.5 Rat Monoclonal Antibody [Clone ID: KT44]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	KT44
Applications:	ELISA, IF, WB
Recommended Dilution:	Western blotting: 1/500-1/5000 (Use 1/2000 as a starting dilution). Immunofluorescence: 1/100-1/1000 (Use 1/200 as a starting dilution).
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	GST fusion protein having the APC region A (788-1038) shown in Figure 1.
Specificity:	Recognizes Mouse Adenomatous Polyposis coli protein (APC).
Formulation:	PBS containing 0.09% Sodium Azide as a preservative State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Protein G Chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for three months or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Database Link:	<a href="#">Q61315</a>
Synonyms:	Adenomatous polyposis coli protein



[View online »](#)

Product images:




Figure 1. Mouse APC