

## Product datasheet for **AM20828PU-L**

### CLASP2 (N-term) Rat Monoclonal Antibody [Clone ID: KT68]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	KT68
Applications:	IF, WB
Recommended Dilution:	Western blotting: 1 µg/ml. Immunofluorescent staining: 10 µg/ml.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	alpha-CLASP2-N-terminus and GST fusion protein
Specificity:	This antibody reacts to CLASP2 (N-term).
Formulation:	PBS containing 0.09% sodium azide as a preservative State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Protein G chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	CLIP associating protein 2
Database Link:	<a href="#">Entrez Gene 76499 Mouse Q8BRT1</a>
Synonyms:	CLIP-associating protein 2, KIAA0627, hOrbit2



[View online »](#)