

Product datasheet for **AM20829PU-N**

Rat IgG Light Chain Mouse Monoclonal Antibody [Clone ID: KT100]

Product data:

| | |
|-----------------------|---|
| Product Type: | Secondary Antibodies |
| Product Name: | Rat IgG Light Chain Mouse Monoclonal Antibody [Clone ID: KT100] |
| Clone Name: | KT100 |
| Applications: | ELISA, WB |
| Recommended Dilution: | ELISA: 1/3000. Western Blot: 1/1000 dilution using DAB as substrate, 1/10000 dilution using ECL system. |
| Reactivity: | Rat |
| Host: | Mouse |
| Immunogen: | Rat IgG |
| Isotype: | IgG2b |
| Formulation: | PBS containing 0.01% thimerosal as a preservative State: Aff - Purified State: Liquid purified Ig fraction Label: Horseradish Peroxidase (HRP) |
| Concentration: | lot specific |
| Purification: | Protein A chromatography |
| Conjugation: | Unconjugated |
| Storage: | Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing. |



[View online »](#)