

Product datasheet for **AM20829PU-S**

Rat IgG Light Chain Mouse Monoclonal Antibody [Clone ID: KT100]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG Light Chain Mouse Monoclonal Antibody [Clone ID: KT100]
Clone Name:	KT100
Applications:	ELISA, WB
Recommended Dilution:	ELISA: 1/3000. Western Blot: 1/1000 dilution using DAB as substrate, 1/10000 dilution using ECL system.
Reactivity:	Rat
Host:	Mouse
Immunogen:	Rat IgG
Isotype:	IgG2b
Formulation:	PBS containing 0.01% thimerosal as a preservative State: Aff - Purified State: Liquid purified Ig fraction Label: Horseradish Peroxidase (HRP)
Concentration:	lot specific
Purification:	Protein A chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C to -80°C for longer. Avoid repeated freezing and thawing.



[View online »](#)