

Product datasheet for **AM20842PU-N**

Prolactin (PRL) Mouse Monoclonal Antibody [Clone ID: n.a]

Product data:

Product Type:	Primary Antibodies
Clone Name:	n.a
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human Prolactin
Specificity:	This antibody reacts to Prolactin.
Formulation:	5 mM Citric Acid pH 2.5, neutralization with 50 mM Na ₂ HPO ₄ at pH 8.8. State: Aff - Purified State: Liquid purified Ig fraction
Purification:	Protein A or Protein G Affinity Chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	prolactin
Database Link:	Entrez Gene 5617 Human P01236
Synonyms:	prolactin
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Cytokine-cytokine receptor interaction, Jak-STAT signaling pathway, Neuroactive ligand-receptor interaction



[View online »](#)