

Product datasheet for **AM20959PU-N**

LMO4 Mouse Monoclonal Antibody [Clone ID: 2B6]


Product data:

Product Type:	Primary Antibodies
Clone Name:	2B6
Applications:	ELISA, IHC
Recommended Dilution:	ELISA. Immunohistochemistry on Paraffin Sections: 5 µg/ml.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Recombinant protein, LMO4 (AAH17673, 1 a.a. ~ 166 a.a) full length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Specificity:	This antibody reacts to LIM Domain Only 4 (LMO4).
Formulation:	PBS, pH 7.2 State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Protein A Chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	LIM domain only 4
Database Link:	Entrez Gene 8543 Human P61968
Synonyms:	LMO4, LIM domain transcription factor LMO4, LIM domain only protein 4, LMO-4, Breast tumor autoantigen
Protein Families:	Transcription Factors



[View online »](#)

Product images:



Human Breast: Formalin-Fixed, Paraffin-Embedded (FFPE)