

## Product datasheet for **AM26014PU-N**

### Human IgE Mouse Monoclonal Antibody [Clone ID: 4G7]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgE Mouse Monoclonal Antibody [Clone ID: 4G7]
Clone Name:	4G7
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Human
Host:	Mouse
Immunogen:	Purified human IgE
Isotype:	IgG
Formulation:	Tris buffered saline (TBS) with 15 mM sodium azide, approx. pH 8.0 State: Purified State: Liquid Ig fraction
Concentration:	lot specific
Purification:	Precipitation methods and ion exchange chromatography (> 95% pure by SDS-PAGE)
Conjugation:	Unconjugated
Storage:	Store at 2 - 8 °C for up to one month or (in aliquots) at -20 °C for longer. Avoid repeated freezing and thawing.



[View online »](#)