

Product datasheet for **AM26023BT-N**

Cd9 Rat Monoclonal Antibody [Clone ID: EM-04]

Product data:

Product Type:	Primary Antibodies
Clone Name:	EM-04
Applications:	FC, WB
Recommended Dilution:	Flow cytometry.
Reactivity:	Mouse
Host:	Rat
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Permeabilized murine bone marrow-derived mast cells (BMMC)
Specificity:	The antibody recognizes CD9 antigen, a 24 kDa transmembrane protein expressed on platelets, monocytes, pre-B lymphocytes, granulocytes and activated T lymphocytes.
Formulation:	Phosphate buffered saline (PBS), approx. pH 7.4 Label: Biotin State: Liquid purified Ig fraction Preservative: 15 mM sodium azide Label: Conjugated with Biotin-LC-NHS under optimum conditions. The reagent is free of unconjugated biotin.
Concentration:	lot specific
Conjugation:	Biotin
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE!
Stability:	Shelf life: one year from despatch.
Gene Name:	CD9 antigen
Database Link:	Entrez Gene 12527 Mouse P40240




[View online »](#)

Background:

CD9 belongs to proteins of tetraspanin family that orchestrate cholesterol-associated tetraspanin-enriched signaling microdomains within the plasma membrane, forming complexes with each other as well as with integrins, membrane-anchored growth factors and other proteins. CD9 is involved in cell motility, osteoclastogenesis, neurite outgrowth, myotube formation, and sperm-egg fusion, plays roles in cell attachment and proliferation and is necessary for association of heterologous MHC II molecules on the dendritic cell plasma membrane which is important for effective T cell stimulation. CD9 is also considered as metastasis suppressor in solid tumors.

Synonyms:

MIC3, TSPAN29, GIG2, p24, Tetraspanin-29, 5H9 antigen

Product images:

Detection of CD9 in bone marrow-derived mast cell (BMMC) lysates by rat monoclonal EM-04 antibody.