

## Product datasheet for **AM26041PU-N**

### Amyloid beta (1-40 specific) Mouse Monoclonal Antibody [Clone ID: CV9 8B2]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	CV9 8B2
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	C-terminal amino acid sequence of human beta amyloid peptide 40, conjugated with KLH
Specificity:	This antibody recognizes the C-terminal peptide of beta amyloid peptide 40, and full length beta amyloid peptide 40.
Formulation:	0.01M PBS pH7.2 State: Aff - Purified State: Lyophilized Ig fraction
Reconstitution Method:	Double distilled water is recommended to adjust the final concentration to 1.0 mg/ml.
Purification:	Protein G affinity chromatography
Conjugation:	Unconjugated
Storage:	Store at (in aliquots) at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)