

Product datasheet for **AM26047SU-N**

Parathyroid Hormone (PTH) Mouse Monoclonal Antibody [Clone ID: MPTH1]

Product data:

Product Type:	Primary Antibodies
Clone Name:	MPTH1
Applications:	ELISA, IHC
Reactivity:	Human
Host:	Mouse
Isotype:	IgM
Clonality:	Monoclonal
Immunogen:	pTH 1-34
Specificity:	This antibody detects Parathyroid Hormone.
Formulation:	Stock solution contains 0.4 mg/ml IgM, phosphate buffered saline pH 7.2 (PBS), 10 mg/ml bovine serum albumin (BSA) as a stabilizer and 0.01 % thimerosal as a preservative. State: Supernatant State: Lyophilized cell culture supernatant
Reconstitution Method:	Reconstitute by adding 0.5ml distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at 2-8°C. Following reconstitution store (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	parathyroid hormone
Database Link:	Entrez Gene 5741 Human P01270
Synonyms:	Parathormone, Parathyrin
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)