

Product datasheet for **AM26086PU-N**

HIV-1 Tat (71-81) Mouse Monoclonal Antibody [Clone ID: N1 (02-012)]

Product data:

Product Type:	Primary Antibodies
Clone Name:	N1 (02-012)
Applications:	ELISA, WB
Recommended Dilution:	Titer with rTatHXB2: ELISA: 1:5000. Western blot: 1:500.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	This antibody detects HIV-1 Tat protein, mapped to amino acids 71-81 (KQPASQPRGDP – B subtype aa. consensus). Reacts with HIV-1 Tat of HXB2, BH10, IIIB and HAN strains.
Formulation:	PBS pH 7.4, with 0.09% sodium azide State: Aff - Purified State: Liquid Ig fraction
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	Unconjugated
Storage:	Upon receipt, store (in aliquots) at -20 to -70 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Background:	HIV is a highly variable virus which mutates very readily. This means there are many different strains of HIV, even within the body of a single infected person. The strains of HIV1 can be classified into three groups: the "major" group M, the "outlier" group O and the "new" group N. These three groups may represent three separate introductions of simian immunodeficiency virus into humans. Group O appears to be restricted to West-Central Africa and group N, discovered in 1998 in Cameroon, is extremely rare. More than 90% of HIV1 infections belong to HIV1 group M.



[View online »](#)

Synonyms: HIV1, HIV-I, Human immunodeficiency virus type 1