

Product datasheet for **AM26140AC-N**

Human Lambda Light Chain Mouse Monoclonal Antibody [Clone ID: 1-155-2]


Product data:

Product Type:	Primary Antibodies
Clone Name:	1-155-2
Applications:	FC
Recommended Dilution:	Flow cytometry: 10 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Human IgA1 lambda myeloma protein
Specificity:	This antibody detects human Lambda light chains.
Formulation:	Phosphate buffered saline (PBS) Label: APC State: Liquid purified Ig fraction Preservative: 15 mM sodium azide
Conjugation:	APC
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE! This products is photosensitive and should be protected from light.
Stability:	Shelf life: one year from despatch.
Background:	Immunoglobulin classes share the same basic four polypeptide chain structure of two heavy chains (five heavy chains types) and two light chains (kappa, lambda; both having a molecular weight of 22.5kDa). Kappa and lambda consist of a variable region and a constant region and can easily be differentiated by the antigenic properties of the constant region. The ratio of kappa to lambda is 70:30.



[View online »](#)

Product images:



Surface staining of lambda light chains in human peripheral blood with anti-lambda light chains (1-155-2) APC.