

Product datasheet for **AM26144FC-N**

CD73 (NT5E) Mouse Monoclonal Antibody [Clone ID: AD2]

Product data:

Product Type:	Primary Antibodies
Clone Name:	AD2
Applications:	FC
Recommended Dilution:	Flow Cytometry analysis of Human blood cells using 4 µl reagent / 100 µl of whole blood or 10 ⁶ cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	This antibody AD2 recognizes CD73, a 70 kDa GPI-anchored 5' -nucleotidase expressed predominantly on T and B cell subsets, follicular dendritic cells and endothelial cells.
Formulation:	PBS Label: FITC State: Liquid purified Ig fraction Preservative: 15 mM Sodium Azide Label: Conjugated with Fluorescein isothiocyanate (FITC) under optimum conditions. The reagent is free of unconjugated FITC and adjusted for direct use.
Conjugation:	FITC
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE! This products is photosensitive and should be protected from light.
Stability:	Shelf life: one year from despatch.
Gene Name:	5'-nucleotidase ecto
Database Link:	Entrez Gene 4907 Human P21589




[View online »](#)

Background:

CD73 (ecto-5'-nucleotidase) is a 70 kDa glycoprotein anchored to the extracellular leaflet of the plasma membrane by GPI. This ecto-enzyme catalyzes dephosphorylation of AMP to adenosine. CD73 is expressed in various types of cells, such as epithelial, muscle, and endothelial cells, neutrophils, lymphocytes and fibroblasts. Inflammatory mediators support CD73 expression and its enzymatic activity, leading to the release of adenosine, which modulates inflammation through adenosine receptors. CD73 is expressed in a variety of lymphomas and leukemias, including ALL and CLL, whereas immunodeficient patients usually express low levels of this protein.

Synonyms:

NT5, NTE, NT5E

Product images:

Surface staining (flow cytometry) of human peripheral blood with anti-human CD73 (AD2) FITC.