

## Product datasheet for **AM26289FC-N**

### TNF alpha (TNF) Mouse Monoclonal Antibody [Clone ID: T1]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	T1
Applications:	ELISA, WB
Recommended Dilution:	Flow cytometry: The typical starting working dilution is 1:50. Western blot. Immunoassay.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	The antibody reacts with free soluble (17 kDa) and membrane (26 kDa) human TNF-alpha. The antibody inhibits the biological activity of both forms. It does not react with receptor bound TNF-alpha. It can be a useful tool to discriminate between the membrane form of TNF expressed on producer cells and the proteolytically cleaved, soluble TNF-alpha bound to its cognate cell membrane receptors (TNF-RI and TNF-RII).
Formulation:	PBS Label: FITC State: Liquid 0.2 µm filtered Ig fraction Stabilizer: 1% bovine serum albumin Preservative: 0.02% sodium azide
Concentration:	lot specific
Purification:	Protein G
Conjugation:	FITC
Storage:	Store at 2 - 8 °C.
Stability:	Shelf life: one year from despatch.
Gene Name:	tumor necrosis factor
Database Link:	<a href="#">Entrez Gene 7124 Human P01375</a>



[View online »](#)

**Synonyms:** TNF, TNF-a, TNFA, TNFSF2, Cachectin