

## Product datasheet for **AM26353PU-N**

### Chicken IgG / Chicken IgY Mouse Monoclonal Antibody [Clone ID: 7C2]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Chicken IgG / Chicken IgY Mouse Monoclonal Antibody [Clone ID: 7C2]
Clone Name:	7C2
Applications:	ELISA, IP, WB
Recommended Dilution:	Immunoassays. Western blot (non-reduced): The typical starting working dilution is 1:10. Also useful for immunopurification.
Reactivity:	Chicken, Duck
Host:	Mouse
Isotype:	IgG1
Formulation:	PBS State: Purified State: Liquid 0.2 µm filtered Ig fraction Stabilizer: 0.1% bovine serum albumin Preservative: 0.02% sodium azide
Concentration:	lot specific
Purification:	Protein G
Conjugation:	Unconjugated
Storage:	Store at 2 - 8 °C.



[View online »](#)