

Product datasheet for AM26704PU-N

Mouse IgG2b Isotype Control

Product data:

Product Type: Isotype Controls

C, FC, ICC/IF, IP, P, WB **Applications:**

Host: Mouse Isotype: IgG2b

Clonality: Monoclonal

Immunogen: KLH-coupled Trinitrophenol

Specificity: The reagent is intended as isotype control for flow cytometry analysis to establish the amount

> of non-specific antibody binding. This Mouse IgG2b (kappa) monoclonal antibody (clone MPC-11) reacts with an epitope irrelevant for a variety of resting, activated, live, and fixed Human,

Mouse, and Rat tissues.

Buffer: PBS, pH~7.4 with 15 mM Sodium Azide

Concentration: 1 mg/ml

Purification: Affinity Chromatography on Protein A

Conjugation: Unconjugated

Storage: Store undiluted at 2-8°C.

DO NOT FREEZE!

Shelf life: one year from despatch.

Background: The specificity of staining by monoclonal antibodies to target antigens should be verified by

> establishing the amount of non-specific antibody binding. Especially at higher concentration (more than 15 µg/ml) the antibody staining usually has consignable background. To this end a non-reactive immunoglobulin of the same isotype is included as a negative control for each specific monoclonal antibody used in a particular immunoassay. The monoclonal antibody MPC-11 was generated against an epitope irrelevant for human, mouse, and rat material, and

can thus be used for evaluation of the background staining that is caused by general nonspecific interactions between an mouse IgG2b molecule and the respective sample under the particular conditions. This shall help the customer to set up the experimental conditions

so that the nonspecific binding of any antibody is abolished.



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com