

## Product datasheet for **AM26775PU-N**

### Mouse IgG2a Isotype Control

#### Product data:

Product Type:	Isotype Controls
Applications:	FC, ICC/IF, WB
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Mineral oils
Specificity:	This antibody reacts with an unknown epitope. It does not react with a variety of resting, activated, live, and fixed mouse, rat and human tissues.
Buffer:	Phosphate buffered saline (PBS)
Concentration:	1 mg/ml
Purification:	Protein-A affinity chromatography (>95% pure by SDS-PAGE)
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C. <b>DO NOT FREEZE!</b> Shelf life: one year from despatch.
Background:	The specificity of staining by monoclonal antibodies to target antigens should be verified by establishing the amount of non-specific antibody binding. Especially at higher concentration (more than 15 µg/ml) the antibody staining usually has consignable background. To this end a non-reactive immunoglobulin of the same isotype is included as a negative control for each specific monoclonal antibody used in a particular immunoassay. The monoclonal antibody MOPC-173, generated against an undefined antigen, does not react specifically with mouse, rat and human samples, and hence all the background that could be observed when working with this antibody would be a result of general nonspecific interactions between an mouse IgG2a molecule and the respective sample under the particular conditions. This shall help the customer to set up the experimental conditions so that the nonspecific binding of any antibody is abolished.



[View online »](#)