

## Product datasheet for **AM26793PU-N**

### Human IgD Mouse Monoclonal Antibody [Clone ID: IA6-2]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgD Mouse Monoclonal Antibody [Clone ID: IA6-2]
Clone Name:	IA6-2
Applications:	FC
Recommended Dilution:	<b>Flow cytometry.</b>
Reactivity:	Human
Host:	Mouse
Immunogen:	Human IgD
Isotype:	IgG2a
Formulation:	Phosphate buffered saline (PBS) State: Purified State: Liquid Ig fraction Preservative: 15 mM sodium azide, approx. pH 7.4
Concentration:	lot specific
Purification:	Protein-A affinity chromatography (> 95% pure by SDS-PAGE)
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C. <b>DO NOT FREEZE!</b>



[View online »](#)