

Product datasheet for AM26793RP-N

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Human IgD Mouse Monoclonal Antibody [Clone ID: IA6-2]

Product data:

Product Type: Secondary Antibodies

Product Name: Human IgD Mouse Monoclonal Antibody [Clone ID: IA6-2]

Clone Name: IA6-2 Applications: FC

Recommended Dilution: Flow cytometry analysis of human blood cells using 10 μl reagent / 100 μl of whole blood or

106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.

Reactivity: Human
Host: Mouse

Immunogen: Human IgD

Isotype: IgG2a

Formulation: Phosphate buffered saline (PBS)

Label: PE

State: Liquid Ig fraction

Stabilizer: 0.2% (w/v) high-grade protease free Bovine Serum Albumin (BSA)

Preservative: 15 mM sodium azide

Label: Conjugated with R-Phycoerythrin (PE) under optimum conditions. The conjugate is

purified by size-exclusion chromatography and adjusted for direct use.

Conjugation: PE

Storage: Atore the antibody undiluted at 2-8°C.

DO NOT FREEZE!

This product is photosensitive and should be protected from light.

