

Product datasheet for **AM28000HR-N**

Goat IgG (Fc specific) Mouse Monoclonal Antibody [Clone ID: SB115d]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Goat IgG (Fc specific) Mouse Monoclonal Antibody [Clone ID: SB115d]
Clone Name:	SB115d
Applications:	ELISA
Recommended Dilution:	ELISA: 1/10,000-1/20,000.
Reactivity:	Goat
Host:	Mouse
Immunogen:	Pooled Goat IgG
Isotype:	IgG2a
Formulation:	50% Glycerol/50% PBS, pH 7.4; without preservatives Label: HRP State: Liquid purified Ig fraction. Label: Horseradish peroxidase
Conjugation:	HRP
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)