

## **Product datasheet for AM31061PU-N**

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## PC1/3 (PCSK1) (652-754) Mouse Monoclonal Antibody [Clone ID: 3D2]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: 3D2

**Applications:** ELISA, IF, IHC, WB

Recommended Dilution: ELISA.

**Immunofluorescence:** 30 µg/ml.

Immunohistochemistry on Paraffin Sections: 10 µg/ml.

Western Blot: 1/500 - 1/1000.

Reactivity: Human, Mouse, Rat

Host: Mouse Isotype: IgG3

Clonality: Monoclonal

**Immunogen:** PCSK1 antibody was raised against against PCSK1 (NP\_000430, 652 a.a. ~ 754 a.a) partial

recombinant protein with GST tag.

**Specificity:** This antibody reognizes PCSK1.

**Formulation:** PBS, pH 7.4

State: Purified

State: Liquid purfied Ig fraction

**Concentration:** lot specific

**Purification:** Protein A Chromatography

Conjugation: Unconjugated

Storage: Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.

**Stability:** Shelf life: one year from despatch.

**Gene Name:** proprotein convertase subtilisin/kexin type 1

Database Link: Entrez Gene 25204 RatEntrez Gene 5122 Human

P29120

**Background:** PCSK1 is a Serine S8 type protease.






Synonyms: PCSK1, Prohormone convertase 1, Prohormone convertase 3, Proprotein convertase 1, NEC1,


NEC 1, Neuroendocrine convertase 1, BMIQ12

**Protein Families:** Druggable Genome, Protease, Secreted Protein


## **Product images:**



PCSK1 monoclonal antibody clone 3D2. Western Blot analysis of PCSK1 expression in PC-12.




PCSK1 monoclonal antibody clone 3D2. Western Blot analysis of PCSK1 expression in NIH/3T3.




Detection limit for recombinant GST tagged PCSK1 is approximately 0.03ng/ml as a capture antibody.





Human Brain, Cortex: Formalin-Fixed, Paraffin-Embedded (FFPE)



Immunofluorescence of monoclonal antibody to PCSK1 on HeLa cell. [antibody concentration 30 ug/ml]