

Product datasheet for **AM31173FC-N**

IL6ra Mouse Monoclonal Antibody [Clone ID: B-N12]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B-N12
Applications:	ELISA, FC
Recommended Dilution:	Flow Cytometry: Use 10µl to label 10e6 cells or 100µl of whole blood
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	gp80 transfected CHO cell line
Specificity:	Recognises the Interleukin-6 Receptor a subunit, IL-6Ra, IL-6R, a 80kDa protein.
Formulation:	Phosphate-buffered saline containing 5% BSA as stabilizer and 0.09% sodium azide as preservative Label: FITC State: Liquid purified Ig fraction
Purification:	Ion exchange chromatography
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	interleukin 6 receptor, alpha
Database Link:	Entrez Gene 16194 Mouse P22272
Synonyms:	Interleukin-6 receptor subunit alpha, IL-6R-alpha, IL-6R 1, IL6 Receptor alpha



[View online »](#)