

Product datasheet for **AM31188AF-N**


CD56 (NCAM1) Mouse Monoclonal Antibody [Clone ID: B-A19]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B-A19
Applications:	FC
Recommended Dilution:	Flow Cytometry.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	U937 cell line
Specificity:	Recognises the neuronal cell adhesion molecule (N-CAM), a 175-220 kDa protein.
Formulation:	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free State: Azide Free State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Ion exchange chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	neural cell adhesion molecule 1
Database Link:	Entrez Gene 4684 Human P13591
Synonyms:	NCAM-1, N-CAM-1, NCAM
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Transmembrane
Protein Pathways:	Cell adhesion molecules (CAMs), Prion diseases



[View online »](#)

Product images:

A typical staining pattern with the B-A19 monoclonal antibody of lymphocytes