

## Product datasheet for **AM31195PU-N**


### CD95 (FAS) Mouse Monoclonal Antibody [Clone ID: B-G27]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	B-G27
Applications:	FC, IHC
Recommended Dilution:	<b>Immunohistochemistry on frozen sections.</b> <b>Flow Cytometry:</b> Use 10 µl to label 10e6 cells or 100 µl of whole blood
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Recombinant Fas antigen
Specificity:	Recognises the Fas, APO-1 antigen, a 45 kDa protein.
Formulation:	Phosphate-buffered saline containing 1 % BSA as stabilizer and 0.09 % sodium azide as preservative State: Purified State: Liquid purified Ig fraction
Purification:	Ion exchange chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	Fas cell surface death receptor
Database Link:	<a href="#">Entrez Gene 8794 Human</a> <a href="#">Entrez Gene 355 Human</a> <a href="#">P25445</a>
Synonyms:	FASLG receptor, Apo-1 antigen, APT1, FAS1, TNFRSF6
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Apoptosis, Cytokine-cytokine receptor interaction, Natural killer cell mediated cytotoxicity



[View online »](#)

**Product images:**

A typical staining pattern with the B-G27 monoclonal antibody of lymphocytes