

Product datasheet for **AM31217RP-N**


CD5 Mouse Monoclonal Antibody [Clone ID: B-B8]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B-B8
Applications:	FC
Recommended Dilution:	Flow Cytometry: Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	PHA activated T cells
Specificity:	This antibody recognises a 67 kDa protein.
Formulation:	Phosphate-buffered saline containing 1% BSA as stabilizer and 0.09% sodium azide as preservative Label: PE State: Lyophilized purified IgG fraction Label: R - Phycoerythrin
Reconstitution Method:	Restore with 1.0 ml deionised water.
Purification:	Ion exchange chromatography
Conjugation:	PE
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE!
Stability:	Shelf life: 6 month from despatch.
Gene Name:	CD5 molecule
Database Link:	Entrez Gene 921 Human P06127
Synonyms:	CD5, LEU1
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Hematopoietic cell lineage



[View online »](#)

Product images:

A typical staining pattern with the B-B8 monoclonal antibody of lymphocytes