

## Product datasheet for **AM31243AF-N**


### CD151 Mouse Monoclonal Antibody [Clone ID: B-E31]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	B-E31
Applications:	FC
Recommended Dilution:	Flow Cytometry.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Eahy 926 cell line
Specificity:	This antibody recognises the antigen PETA-3, a 32 kDa protein.
Formulation:	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free State: Azide Free State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Ion exchange chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	CD151 molecule (Raph blood group)
Database Link:	<a href="#">Entrez Gene 977 Human P48509</a>
Synonyms:	PETA-3, PETA3, GP27, Tetraspanin-24, SFA-1, SFA1
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

**Product images:**

A typical staining pattern with the B-E31 monoclonal antibody of Eahy 926 cell line