

Product datasheet for **AM31244AF-N**


CD109 Mouse Monoclonal Antibody [Clone ID: B-E47]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B-E47
Applications:	FC
Recommended Dilution:	Flow Cytometry.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Eahy 926 cell line
Specificity:	This antibody recognises a 170 kDa protein.
Formulation:	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free State: Azide Free State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Ion exchange chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	CD109 molecule
Database Link:	Entrez Gene 135228 Human Q6YHK3
Synonyms:	p180, r150



[View online »](#)

Product images:

A typical staining pattern with the B-E47 monoclonal antibody of Eahy 926 cell line