

Product datasheet for **AM31246AF-N**


LFA3 (CD58) Mouse Monoclonal Antibody [Clone ID: B-L28]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B-L28
Applications:	FC
Recommended Dilution:	Flow Cytometry.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	CD58 transfected CHO cell line
Specificity:	This antibody recognises the Lymphocyte Function-associated Antigen-3 (LFA-3), a 40-70 kDa protein.
Formulation:	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free State: Azide Free State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Ion exchange chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	CD58 molecule
Database Link:	Entrez Gene 965 Human P19256
Synonyms:	Ag3, LFA-3
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Cell adhesion molecules (CAMs)



[View online »](#)

Product images:

A typical staining pattern with the B-L28 monoclonal antibody of lymphocytes