

## Product datasheet for **AM31246PU-N**


### LFA3 (CD58) Mouse Monoclonal Antibody [Clone ID: B-L28]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	B-L28
Applications:	FC
Recommended Dilution:	Flow Cytometry: Use 10 µl to label 10e6 cells or 100 µl of whole blood.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	CD58 transfected CHO cell line
Specificity:	This antibody recognises the Lymphocyte Function-associated Antigen-3 (LFA-3), a 40-70 kDa protein.
Formulation:	Phosphate-buffered saline containing 1 % BSA as stabilizer and 0.09 % sodium azide as preservative State: Purified State: Liquid purified Ig fraction
Purification:	Ion exchange chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	CD58 molecule
Database Link:	<a href="#">Entrez Gene 965 Human P19256</a>
Synonyms:	Ag3, LFA-3
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Cell adhesion molecules (CAMs)



[View online »](#)

**Product images:**

A typical staining pattern with the B-L28 monoclonal antibody of lymphocytes