

Product datasheet for **AM31254AF-N**


HLA Class II DQ Mouse Monoclonal Antibody [Clone ID: B-K27]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B-K27
Applications:	FC
Recommended Dilution:	Flow Cytometry.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	T4+ CTL clone HG-38
Specificity:	This antibody recognises HLA Class II DQ a 28kDa and 33kDa complex.
Formulation:	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free State: Azide Free State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Ion exchange chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Synonyms:	MHC Class II DQ



[View online »](#)

Product images:

A typical staining pattern with the B-K27 monoclonal antibody of lymphocytes