

## Product datasheet for **AM31257FC-N**

### Interferon gamma (IFNG) Mouse Monoclonal Antibody [Clone ID: B-B1]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	B-B1
Applications:	FC
Recommended Dilution:	<b>Intra-cytoplasmic Flow cytometry:</b> Use 10 µl of antibody to label 5x10e5 cells.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Recombinant human IFN-gamma
Specificity:	This antibody recognises both natural and recombinant IFN-gamma.
Formulation:	Phosphate-buffered saline containing 5 % BSA as stabilizer and 0.09 % sodium azide as preservative Label: FITC State: Liquid purified Ig fraction Label: Fluoresceinisothiocyanat
Purification:	Ion exchange chromatography
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	interferon, gamma
Database Link:	<a href="#">Entrez Gene 3458 Human P01579</a>
Synonyms:	IFN-gamma, gamma IFN
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)

**Protein Pathways:**

Allograft rejection, Cytokine-cytokine receptor interaction, Graft-versus-host disease, Jak-STAT signaling pathway, Natural killer cell mediated cytotoxicity, Proteasome, Regulation of autophagy, Systemic lupus erythematosus, T cell receptor signaling pathway, TGF-beta signaling pathway, Type I diabetes mellitus