

Product datasheet for **AM31268AF-N**

IL15 Mouse Monoclonal Antibody [Clone ID: B-E29]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B-E29
Applications:	FN
Recommended Dilution:	Functional Assay: Blocks IL-15 induced proliferation on Kit 225 cell line.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Recombinant human IL-15
Specificity:	This antibody recognises both natural and recombinant IL-15.
Formulation:	Phosphate-buffered saline. Sterile-filtered through 0.22 µm and treated to remove endotoxins. Carrier and preservative free State: Azide Free State: Liquid purified IgG fraction
Concentration:	lot specific
Purification:	Ion Exchange Chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	interleukin 15
Database Link:	Entrez Gene 3600 Human P40933
Synonyms:	IL-15
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Cytokine-cytokine receptor interaction, Jak-STAT signaling pathway



[View online »](#)