

Product datasheet for **AM31270AF-N**

IL17F Mouse Monoclonal Antibody [Clone ID: B-F60]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B-F60
Applications:	FN
Recommended Dilution:	Functional Assay: Neutralises the bioactivity of IL-17F on HT1080 cell line (measurement of IL-6 production).
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Recombinant human IL-17F
Specificity:	This antibody recognizes both recombinant and natural Human IL-17F.
Formulation:	Phosphate-buffered saline. Sterile-filtered through 0.22 µm and treated to remove endotoxins. Carrier and preservative free State: Azide Free State: Liquid purified IgG fraction
Concentration:	lot specific
Purification:	Ion exchange chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	interleukin 17F
Database Link:	Entrez Gene 112744 Human Q96PD4
Synonyms:	IL24, Cytokine ML-1, Interleukin-24, IL-24
Note:	The clone B-F60 is also available biotinylated as Cat.No. AM31346BT-N.
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)