

Product datasheet for **AM31278AF-N**

Azurocidin (AZU1) Mouse Monoclonal Antibody [Clone ID: B-P37]

Product data:

Product Type:	Primary Antibodies
Clone Name:	B-P37
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	Recombinant human Azurocidin/Cap37/HBP
Specificity:	This antibody recognises both natural and recombinant human HBP.
Formulation:	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free State: Azide Free State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Ion exchange chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	azurocidin 1
Database Link:	Entrez Gene 566 Human P20160
Synonyms:	CAP37, Heparin-binding protein, HBP
Protein Families:	Druggable Genome, Protease



[View online »](#)