

## Product datasheet for **AM31282AF-N**

### IL33 Mouse Monoclonal Antibody [Clone ID: B-L33]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	B-L33
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Recombinant human IL-33/NF-HEV
Specificity:	This antibody recognises both natural and recombinant IL-33.
Formulation:	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free State: Azide Free State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Ion exchange chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Gene Name:	interleukin 33
Database Link:	<a href="#">Entrez Gene 90865 Human O95760</a>
Synonyms:	IL-33, IL1F11, IL-1F11, NFHEV, NF-HEV, C9orf26
Protein Families:	Secreted Protein
Protein Pathways:	Cytosolic DNA-sensing pathway



[View online »](#)